PHARMACEUTICAL CAPSULES

ACGcaps™ GT





For enrobing caplets.

Using ACGcaps™ GT unleashes great opportunities to differentiate your brand in an increasingly crowded OTC market.

Comprising two identically shaped gelatin shells, ACGcaps™ GT allows you to completely enclose caplets (capsule-shaped tablets). A gap can also be left between shells, with the caplet in the very centre of the capsule.

Key features

- / Creates brand differentiation
- / Masks unpleasant tastes and odours of active ingredients and / or excipients
- / Provides a tamperproof shell for caplets
- / Excellent option for product line extensions
- / ACGcaps™ GT Plus variant offers full capsule enrobing
- / ACGcaps™ GT Pro variant leaves a gap between shells for more rapid disintegration and quick release of drugs
- / Smooth capsule surfaces enhance swallowability



Technical specifications

- / Moisture content: 15.5% to 17.0%
- / Disintegration time in water at 37 \pm 2°C: maximum 15 minutes (with guided disc)



This certification is applicable to certain colours and $\mbox{\sc /}$ or variants.

